


Instructions for Use

i - DIGITS

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	15
FR	Notice d'utilisation	27
ES	Instrucciones para el uso	39
IT	Istruzioni per l'uso	51
NO	Bruksanvisning	63
DA	Brugsanvisning	74
SV	Bruksanvisning	85
FI	Käyttöohjeet	96
EL	Οδηγίες Χρήσης	107
NL	Gebruiksaanwijzing	119
PT	Instruções de Utilização	131
PL	Instrukcja użytkowania	143
CS	Návod k použití	155
TR	Kullanım Talimatları	167
RU	Инструкция по использованию	178
JA	取扱説明書	190
ZH	中文说明书	200
KO	사용 설명서	210
AR	تعليمات الاستخدام	221

ENGLISH

INSTRUCTIONS FOR USE

i-DIGITS

The i-Digits is referred to as “*the device*” in the following document. This document provides information on the indications for use, and handling of the device. It is intended for the user of the device. The device may only be fitted and configured by a qualified practitioner authorized by Össur after completing the corresponding training.

These “*Instructions for Use*” relate to: i-Digits Quantum, i-Digits Access, i-Digits France, analog electrodes and mains charger.

PRODUCT DESCRIPTION

i-Digits is a fully customized partial hand prosthesis with individually powered digits. i-Digits work in conjunction with remaining fingers (**Fig. 1a**). The wristband houses the control electronics and the removable batteries. It is fastened using a buckle and is designed to fit loosely around the distal forearm. The product label can be found on the wristband (**Fig. 1b**).

i-Digits Quantum has 20 different grip options available and 12 customizable my grips. All other i-Digits devices have 12 grip options available. The grips can be accessed and programmed through the My i-Limb app. The available grips for the device are displayed on the quick grips page of the My i-Limb app.

Features comparison			
Control options	i-Digits Quantum	i-Digits France	i-Digits Access
App control	Yes	Yes	Yes
Muscle control	Yes	Yes	Yes
Proximity control	Yes	Yes	-
Gesture control	Yes	-	-
Grips available	20	12	12
my grips	12	-	-

i-Digits is fitted to a socket interface assembled by the prosthetist. Minor redness may be experienced on the residual hand but should disappear shortly after removing the device.

Contact your clinician if the socket is causing discomfort or preventing proper use of the i-Digits.

i-Digits should be fitted by a certified prosthetist.

INDICATIONS FOR USE

- Upper limb amputation
- Congenital upper limb absence

CONTRAINDICATIONS

None known.

INTENDED USE

The i-Digits are intended as a part of a prosthetic system that replaces function of a missing upper limb.

REQUIRED DEVICES

The My i-Limb and Biosim apps can be downloaded from the Apple Store. The two apps require an Apple iOS device supported by the manufacturer, e.g. an iPhone or iPad. See the Apple Store for device compatibility.

SAFETY INSTRUCTIONS

Warnings

i-Digits:

- The end user is the intended operator of the device and is responsible for its use.
- Improper handling or adjustment of the i-Digits may cause malfunction of the device.
- The i-Digits do not provide sensation, heat and moisture cannot be felt. The i-Digits is for low to moderate impact activities only.
- Do not use without approved cover(s).
- Do not use with damaged cover(s).
- Do not disassemble componentry or modify in any way.
- Do not service or perform maintenance when in use.
- Do not carry objects using only the tips of the digits. Carry objects by evenly distributing weight across the digits, close to the knuckles and palm of the hand (**Fig. 2**).
- Do not use for heavy lifting.
- Do not use with machinery with moving parts that may cause personal injury or damage.
- Do not use for extreme activities that may cause injury to a natural hand.
- Do not expose to vibration.
- Do not expose to excessive or high forces, particularly on the fingertips and on the side of the digits.
- Do not expose to water.
- Do not expose to excessive moisture, liquids, dust, high temperatures, or shock.
- Do not use in hazardous environments.
- Do not expose to flames.
- Do not use in or expose to explosive atmospheres.
- The electrode is an APPLIED PART.
- This device is for single patient use.
- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.
- The electrodes may contain nickel.

Batteries:

- Do not bend or exert excessive pressure on the battery.
- Do not pierce the battery.
- Do not disassemble battery.
- Do not expose batteries to high temperatures.
- Do not incinerate batteries.
- Do not short circuit the battery.
- Do not store batteries inside a vehicle.
- Dispose of batteries in accordance with US, European or local regulations.

Precautions

i-Digits:

- Users must comply with local regulations on the operation of automobiles, aircraft, sailing vessels of any kind and any other motorized vehicle or device. It is entirely the user's responsibility to seek confirmation that they are physically and legally able to drive using the i-Digits and to the fullest extent permitted by law.
- Only use with approved Össur accessories and tooling.
- Maintenance, repairs and upgrades may only be performed by qualified Össur technicians and technical partners. Össur will provide upon request information to assist service personnel in repair of device.
- Do not use an i-Digits device to operate electronic devices connected to a mains outlet, as this may impact function.
- It is not recommended to use your device while in close proximity of other medical electrical equipment.

- Do not use while charging is in progress.
- Only use with covers supplied by Össur.
- Always use with covers to avoid risk of electrostatic build-up and discharge.
- Do not use oil-based lotions on the skin, e.g. Vaseline.
- Do not expose electrode to dirt or fluids.

Batteries:

- Only use Össur batteries with this device.
- Only use the Össur charger to charge Össur batteries.
- Make sure the battery is not subject to continued pressure once fitted.
- Batteries are to be replaced annually.
- If the battery has visibly ballooned or swelled:
 - discontinue the charging process immediately
 - disconnect the battery
 - remove to a safe area
 - leave and observe for 15 minutes
 - replace with new battery
 - do not re-use
 - dispose of any leaking batteries in an appropriate manner
- If the device will not be used for a long period of time, it is advisable to remove the battery from the prosthesis.

POWER

Batteries

The device can be employed for up to 600 full opening and closing cycles from a single charge of a set of 800mAh batteries. The number of open and close cycles achieved will vary depending on battery age and device use. After one set is drained, the batteries can be easily swapped to the back-up set.

Batteries are placed into the battery case on the wristband, where they are secured by the holding catch located towards the wrist. Release the holding catch to remove the battery (**Fig. 3**).

Each battery is equipped with an LED indicating the charge level of the battery. When the battery charge level is low, the LED will light up red. The LED will stay lit until the battery has been charged to an adequate level.

Charging the Batteries

Charging time: 90 minutes to 3 hours.

Remove the batteries from the wristband and insert into the charger base unit. Connect the charger base unit with the power cable. Plug the power cable into the power outlet.

The charging state is indicated on the back of the base unit (**Fig. 4**):

- Middle light on: Charger is plugged in.
- 2nd and 5th lights flashing green: batteries are charging.
- 2nd and 5th lights solid green: batteries are charged.
- 1st and 4th red lights on: battery fault, unplug and try again. If lights continue to illuminate, contact your clinician.

Turning On/ Off

- Turn the device on by pressing the wristband button; the LED will light up for several seconds to indicate power is on.
- To turn off, press the button again; the LED will light up for a short moment.

CONTROLLING THE I-DIGITS

There are various control options for accessing automated grips on i-Digits devices. Control options vary between different models.

Identifying the Device Number:

Each i-Digits has a unique device number for identifying the device located on the wristband underneath the battery (Fig. 5).

When connecting your i-Digits device to the My i-Limb app, your device number will be displayed on the connection screen. Selecting the number connects the app to your i-Digits device. Alternatively, when connected to the app your device number is displayed in the 'about' section.

Gesture Control (available only on i-Digits™ Quantum)

Enables an automated grip to be accessed through a smooth motion of the prosthesis in one of four directions (forward, back, left or right). You can customize the grips programmed for each direction using the My i-Limb app.

To access gesture control:

1. Device must be in full open hand mode
2. Hold arm parallel to the ground (elbow bent to 90°)
3. Maintain an open signal until finger twitches
4. Move the hand within 1 second in direction set to access a desired grip
5. i-Digits Quantum will adopt the grip

NOTE: Gesture control can also be accessed using co-contraction, depending on the set-up of the i-Digits through the prosthetist.

App Control

An automated grip can be accessed by touching an icon in the My i-Limb app. These automated grips are called quick grips. i-Digits will exit the grip when the icon is tapped again or by selecting another grip icon to enter.

Muscle Control

Triggers are specific muscle signals that can be used to access an automated grip. Four triggers are available: hold open, double impulse, triple impulse and co-contraction.

You can use the app to activate and program muscle control.

Proximity Control (not available on i-Digits Access)

Grip Chips are small Bluetooth devices which activate grips on the prosthesis when the prosthesis is moved close to them or by tapping the Grip Chip.

Before operating a Grip Chip for the first time, pull the plastic tab out from the device to engage the battery.

To use Grip Chips, make sure the i-Digits device is not connected to My i-Limb app.

- **Proximity** is accessed by fully opening the i-Digits near (15 cm/ 6") to a Grip Chip and then relaxing the muscles. Wait until the grip is accessed, which may take up to 3 seconds. To exit a grip accessed via proximity control, give a long open signal.
- **Double tapping** a Grip Chip activates the grip. Quickly tap a Grip Chip twice, just like double-clicking a computer mouse. The LED on the Grip Chip flashes once when tapping is successful.

NOTE: a 3 second pause between each double tap is required. The pause prevents the grip chip from incorrectly detecting multiple taps in a very short period of time. This can cause the hand to enter a grip and then immediately exit it again.

The i-Digits must be fully open/digits stalled open for proximity or tapping to be successful.

The individual Grip Chips are programmed using the My i-Limb app and can be reprogrammed by the user at any time.

Replacing the Grip Chip battery

The battery is expected to last 6 months. Battery life is reduced if Grip Chip is stored close to the i Digits device.

Battery is not rechargeable.

- To replace the battery, insert fingernail between upper and lower caps and run around joint to release.
- Remove battery. Insert new battery under retaining clip, ensuring writing on battery is facing outwards.
- Replace the caps by pressing upper and lower caps together.

Battery type: CR1616. Replacement batteries are available from customer support upon request.

CLEANING

The i-Digits wristband, digits, cover and electrode surface can be cleaned with a soft damp cloth and mild soap. Isopropyl alcohol can also be used to help with disinfection.

Clean the electrode surface regularly.

Do not use any strong chemicals.

NOTE: Do not submerge the i-Digits, wristband, cover, electrode surface or charger in water.

MAINTENANCE

- Always turn off the i-Digits when not in use.
- Charge the battery after use.
- Ensure the i-Digits is serviced every 12 months.
- For all service or maintenance, contact your prosthetist.

TROUBLESHOOTING

Contact your clinician should any problems occur affecting the use of the i-Digits.

TECHNICAL SPECIFICATIONS

i-Digits	
Voltage	7.4 V (nominal)
Max. current	5 A
Battery capacity	Rechargeable lithium polymer 7.4 V (nominal); 800 mAh capacity
Max device load	20 kg/ 44lbs
Finger carry load (static limit)	5 kg/ 11lbs
Time from full open to full close	0.8 seconds
Expected service life	5 years

Power Supply	
Manufacturer	XP Power
Model No.	VEP36US12
Input	100-240Vac, 50-60Hz, 0.9A Max
Output	12VDC, 3A

ENVIRONMENTAL CONDITIONS

Do not use, transport or store the i-Digits outside of the boundaries in the table below:

	Use	Shipping	Extended Storage
Temperature	0 °C to +40 °C	-40 °C to +70 °C	-25 °C to +70 °C
Relative Humidity	10 % to 100 %	10 % to 100 %	10 % to 100 %
Atmospheric pressure	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

COMPLIANCE

The device complies with:

EN 60601-1/AC

Protection against electrical shock – Class II ME Equipment

Degree of protection against electrical shock – Type BF provides additional protection against electric shock

Degree of protection against ingress of water – IP22 (**IEC 60529**)

EMI/EMC

Complies with standard BS EN 60601-1-2:2007/IEC 60601-1-2

Radio Spectrum Matters (ERM)/Bluetooth

Complies with standard EN 301 489-1 V1.9.2

BLUETOOTH MODULES REGULATED INFORMATION

This device contains the following radio frequency transmitters:

Model	Re	Type and Frequency Characteristics	Effective Radiated Power
Bluetooth Low Energy Dual Mode Module Model BR-LE4.0-D2A	FCC Contains FCC ID: XDULE40-D2 Canada Contains IC: 8456A-LE4D2 Japan Contains transmitter with certificate number  R 205-160268	(Dual Mode) Version V2.1 +ED (GFSK + $\pi/4$ DQPSK + 8DPSK) 2402 -2480 MHz Version V4.0 (GFSK) 2402 -2480 MHz	Adjustable Power (-23 dBm to 10.5dBm) short to long range

ELECTROMAGNETIC COMPATIBILITY

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation. In order to regulate the requirements for Electromagnetic Compatibility (EMC) with the aim to prevent unsafe product situations, the BS EN 60601-1-2/ IEC 60601-1-2 standard has been implemented for all Össur Myoelectric Prosthetic Devices. This standard defines levels of electromagnetic emissions for medical devices.

The i-Digits are suitable for use in any environment except where immersion in water or any other fluid is possible, or where exposure to highly electrical and/or magnetic fields can occur (e.g. electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

Refer to further guidance below regarding the EMC environment in which the device should be used:

Guidance and manufacturer's declaration – electromagnetic emissions		
Össur Myoelectric Prosthetic Devices is intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that it each are used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Not applicable Battery Powered	Össur Myoelectric Prosthetic Devices use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Össur Myoelectric Prosthetic Devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable Battery Powered	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable Battery Powered	

Guidance and manufacturer's declaration – electromagnetic immunity

Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. Portable and mobile RF communications equipment should be used no closer to any part of the Össur Myoelectric Prosthetic Devices , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered No Cables >3m
Surge IEC 61000-4-5	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered
Voltage dips, short interruptions and voltage variations on power supply IEC 61000-4-11	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered No Cables >3m
Power frequency (50/ 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m 50/60Hz	Battery Powered

Guidance and manufacturer's declaration – electromagnetic immunity

Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customers or the users of Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	Not applicable	Not applicable Battery Powered No Cables >3m	Portable and mobile RF communications equipment should be used no closer to any part of the Össur Myoelectric Prosthetic Devices including cables, than the recommended separation distance calculated from the equation appropriate to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	12 V/m 26 MHz to 1000 MHz	12 V/m 26 MHz to 1000 MHz	Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz
	10 V/m 1000 MHz to 2700 MHz	10 V/m 1000 MHz to 2700 MHz 1kHz 80% AM	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b Interference may occur in the vicinity of equipment marked with the following symbol:
			

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Össur Myoelectric Prosthetic Devices are used exceeds the applicable RF compliance level above, Össur Myoelectric Prosthetic Devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Össur Myoelectric Prosthetic Devices

Recommended separation distance between portable and mobile RF communications equipment and the Össur Myoelectric Prosthetic Devices

Össur Myoelectric Prosthetic Devices are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customers or the user of Össur Myoelectric Prosthetic Devices can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Össur Myoelectric Prosthetic Devices as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in Watt	Separation distance according to frequency of transmitter in meters		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

DESCRIPTION OF SYMBOLS

BF Applied Part

Manufacturer - YYYY

Manufacturer and year of manufacturing (YYYY)

Consult Instructions for use

Device contains electronic components and/or batteries that should not be disposed of in regular waste

IP22

Protected against solid foreign objects of 12.5mm diameter and greater and protected against splashing water.

SN

Serial Number

For i-Digits™ devices:

The unique serial number for i-Digits™ devices is 2 letters followed by 6 digits.

European Conformity

MD

Medical Device

Class II equipment – provides double isolation to protect against electric shock

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