
## Instructions for Use

---



i - LIMB® HAND

3

EN   Instructions for Use	5
DE   Gebrauchsanweisung	17
FR   Notice d'utilisation	30
ES   Instrucciones para el uso	43
IT   Istruzioni per l'uso	56
NO   Bruksanvisning	70
DA   Brugsanvisning	82
SV   Bruksanvisning	94
EL   Οδηγίες Χρήσης	106
FI   Käyttöohjeet	119
NL   Gebruiksaanwijzing	131
PT   Instruções de Utilização	144
PL   Instrukcja użytkowania	157
CS   Návod k použití	170
TR   Kullanım Talimatları	183
RU   Инструкция по использованию	195
JA   取扱説明書	208
ZH   中文说明书	220
KO   사용 설명서	231
AR   تعليمات الاستخدام	243

1a


1b


LUMI - Orsi HD Grippers 1, 3, 5, 10-Brackets, Inclined PTFE, IP67  
Contains FCC ID 2ABAL-LD022 & IC:2494A-LD022 Made in UK  
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.


2

3


4

5

6


7

8

9


10


11



12

13

**14****16****15****17****18****19****20****21****22**

# ENGLISH

---

## INSTRUCTIONS FOR USE

---

### I-LIMB HAND

The i-Limb® hand is referred to as “*the device*” in the following document. This document provides information on the indications for use, and handling of the device. It is intended for the user of the device. The device may only be configured and fitted by a qualified practitioner authorized by Össur after completing the corresponding training.

These “*Instructions for Use*” relate to: i-Limb Quantum, i-Limb Ultra, i-Limb Access, analog electrodes, magnetic charger port, mains and car charger.

## PRODUCT DESCRIPTION

---

The i-Limb hands are a range of prosthetic hands consisting of individual motorized digits, stall detection and a proprietary control app for Apple iOS devices (**Fig. 1a**). The product label can be found on the side of the device (**Fig. 1b**). The device is assembled as part of your prosthesis together with a socket custom made for you by your clinician.

Users can choose from a selection of automated grips and gestures to help complete their daily tasks. Depending on model, grips can be customized and automated for additional control. An overview of available control options can be found in the features comparison table below:

Features comparison			
Control options	i-Limb Quantum	i-Limb Ultra	i-Limb Access
Gesture control	Yes	-	-
App control	Yes	Yes	Yes
Muscle control	Yes	Yes	Yes
Proximity control	Yes	-	-
Grips available	24	18	12
My grips	12	-	-
Speed boost	Yes	Yes	Yes
Powered thumb rotation	Yes	Yes	-
vari-grip	Yes	Yes	-
Natural hand mode	Yes	Yes	-

While digits are not controlled individually, automated grips allow the i-Limb hand user to move certain digits for pinching an object or making a gesture. Training is required to fully utilize the i-Limb hand and fully understand all the functional benefits.

## INDICATIONS FOR USE

---

- Upper limb amputation
- Congenital upper limb absence

## CONTRAINDICATIONS

---

None known.

## INTENDED USE

---

The i-Limb hand is intended as a part of a prosthetic system that replaces function of a missing upper limb.

## REQUIRED DEVICES

---

The My i-Limb and Biosim apps can be downloaded from the Apple Store. The two apps require an Apple iOS device supported by the manufacturer, e.g. an iPhone or iPad. See the Apple Store for device compatibility.

## SAFETY INSTRUCTIONS

---

### **Warnings**

#### **i-Limb Hand:**

- The end user is the intended operator of the device and is responsible for its use.
- Improper handling or adjustment of the i-Limb hand may cause malfunction of the device.
- The i-Limb hand does not provide sensation, heat and moisture cannot be felt. The i-Limb hand is for low to moderate impact activities only.
- Do not use without an approved cover.
- Do not use with a damaged cover.
- Do not disassemble componentry or modify in any way.
- Do not service or perform maintenance when in use.
- Do not carry objects using only the tips of the digits. Carry objects by evenly distributing weight across the digits, close to the knuckles and palm of the hand.
- Do not use for heavy lifting.
- Do not use with machinery with moving parts that may cause personal injury or damage.
- Do not use for extreme activities that may cause injury to a natural hand.
- Do not expose to vibration.
- Do not expose to excessive or high forces, particularly on the fingertips and on the side of the digits.
- Do not expose to water.
- Do not expose to excessive moisture, liquids, dust, high temperatures, or shock.
- Do not use in hazardous environments.
- Do not expose to flames.
- Do not use in or expose to explosive atmospheres.
- The electrode(s), magnetic charger port, DC port and switch block are an APPLIED PART.
- This device is for single patient use.
- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.
- The electrode may contain nickel.

#### **Batteries:**

- Do not cut or modify battery wires.
- Do not bend or exert excessive pressure on the battery.
- Do not pierce the battery.
- Do not disassemble battery.
- Do not expose batteries to high temperatures.
- Do not incinerate batteries.
- Do not alter battery terminal wires.
- Do not short circuit the battery.
- Do not store batteries inside a vehicle.
- Dispose of batteries in accordance with US, European or local regulations.

### **Precautions**

#### **i-Limb Hand:**

- The user must comply with local regulations on the operation of automobiles, aircraft, sailing vessels of any kind and any other motorized vehicle or device. It is entirely the user's responsibility to seek confirmation that they are physically and legally able to drive using the i-Limb hand and to the fullest extent permitted by law.
- Only use with approved Össur accessories and tooling.
- Damaged covers must be replaced or repaired by a trained individual. Maintenance, repairs and upgrades may only be performed by qualified Össur technicians and technical partners. Össur will provide upon request information to assist service personnel in repair of device.

- Using an i-Limb device to operate electronic devices connected to a mains outlet may impact function.
- It is not recommended to use your device while in close proximity of other medical electrical equipment.
- Do not use while charging is in progress.
- Only use with gloves supplied by Össur.
- Always use with gloves to avoid risk of electrostatic build-up and discharge.
- Do not use oil-based lotions on the skin, e.g. Vaseline.
- Do not expose electrode to dirt or fluids.

#### **Batteries:**

- Only Össur batteries should be used with this device.
- Only use the Össur charger to charge Össur batteries.
- Internal batteries must not be replaced by the end user.
- Batteries are to be replaced annually by service personnel only.
- Do not use your i-Limb hand if the battery has visibly ballooned or swelled. Contact your clinician and do as follows:
  - discontinue the charging process immediately
  - disconnect the battery
  - remove to a safe area
  - leave and observe for 15 minutes
  - replace with new battery
  - do not re-use
  - dispose of any leaking batteries in an appropriate manner
- If the device will not be used for a long period of time, it is advisable to remove the battery from the prosthesis. Contact your clinician for this to be completed.

---

## SAFE USE GUIDELINES

---

- Move the thumb to the side so that it doesn't block the index from closing. Hold handles or similarly shaped objects securely in palm near the base of the fingers (**Fig. 2**).
- Do not block the index finger from closing around thinner objects. The object will be less secure within the i-Limb hand if the digits are not able to conform around its shape (**Fig. 3**).
- Hold objects close to palm with all fingers fully closed. Ensure all digits are fully closed around objects (**Fig. 4**).
- Do not hold objects with fingertips or side of the fingers (**Fig. 5**).
- Position objects close to the knuckles when pulling/pushing (**Fig. 6**).
- Do not pull or push objects with fingertips (**Fig. 7**).
- Push up with fully closed fist with force at knuckles (**Fig. 8**).
- Do not push up on fingers (**Fig. 9**).

---

## POWER

---

The device can be employed for up to 1300 full opening and closing cycles from a single charge of a 1300 mAh battery pack. The number of open and close cycles achieved will vary depending on battery age, use with other associated components (i. e. i-Limb Wrist), and device use as tasks applying more load to the device may drain the batteries faster.

To charge the i-Limb hand, remove the prosthesis from the residual limb. Depending on model, turn either off or to standby by moving power switch to the left position (**Fig. 10, 11**).

#### **Magnetic Charger Port**

The magnetic charger port allows the battery to be charged, the power switched on/off and the remaining battery level to be monitored in one place (**Fig. 12**).

If your i-Limb has an internal battery and magnetic charger port, please follow these steps:

#### **Turning the Device On and Off**

- Press switch on charge port for 1 second and release to turn on/off.
- When turned on: the display will fully light up and then dim to a low level.

- When turned off: the display will fully light up and then switch off.

### **Checking the Battery Charge Level Status**

- The battery charge level is shown by the bars when the power is on.
- Each bar of light shows a 20% charge level. The number of bars lit up shows how much power is left in the battery.
- When the battery level reaches 5%, the red low battery warning will light up (**Fig. 13**). The light will shine for 3 minutes, and then the device will turn itself off.

**NOTE:** an emergency reserve of power remains available, to turn on and release the device from an object, ensuring the safety of the user.

### **Charging the Battery**

**WARNING:** Do not wear an i-Limb hand when the battery is charging.

- Plug the mains charger into the power outlet. Attach the magnetic charger to the charger port on the prosthesis.
- When in standby, the charger unit will display a faint green light.
- When charging is in progress, a red light is displayed.
- When fully charged a green light is displayed.
- Charging time: 90 minutes to 3 hours.
- A car charger is also available.

### **External Battery**

If the i-Limb has external batteries, remove the batteries from the prosthesis and insert into the charger base unit (**Fig. 14**). Connect the charger base unit with the power cable. Plug the power cable into the power outlet.

The charging state is indicated on the back of the base unit:

- Middle light on: charger is plugged in
- 2nd and 5th lights flashing green: batteries are charging
- 2nd and 5th lights solid green: batteries are charged
- 1st and 4th red lights on: battery fault, unplug and try again. If lights continue to illuminate, contact your clinician.

### **Power Supply**

**CAUTION:** Only use the Össur power supply to charge Össur batteries.

<b>Manufacturer:</b>	PowerSolve
<b>Model No.:</b>	000311A
<b>Input:</b>	100-240Vac, 50-60 Hz, 0.3A Max
<b>Output:</b>	8.4VDC, 1A

---

## CONTROLLING THE I-LIMB HAND

---

### **Identifying the Device Number:**

The device number can be found proximal to the thumb base (**Fig. 15**). For extra small hands, or those with a flexion wrist, the number is located within the chassis.

When connecting your i-Limb hand with the My i-Limb app, the device number will be displayed on the connection screen. Selecting the number connects the app to your i-Limb hand. Alternatively, when connected to the app your device number is displayed in the ‘about’ section.

## **Gesture Control (available only on i-Limb Quantum hands)**

Enables an automated grip to be accessed through a smooth motion of the prosthesis in one of four directions (forward, back, left or right). The grips programmed for each direction are customized to the user's requirements using the My i-Limb app.

To access gesture control:

1. Hold i-Limb hand parallel to the ground (elbow bent to 90°)
2. Maintain an open signal until index finger twitches
3. Move the i-Limb Quantum hand (within 1 second) in direction assigned to the desired grip
4. The i-Limb Quantum hand will adopt the grip
5. To exit the grip, hold an open signal until the hand exits the grip

Default setting for gesture control is hold open, however it can alternatively be accessed using co-contraction. This is changed using the My i-Limb app.

## **App Control (available on all i-Limb hands)**

An automated grip can be accessed at the touch of an icon on the My i-Limb app. These are called quick grips. The i-Limb hand will exit the grip when the icon is tapped again or by selecting another grip icon to enter.

## **Muscle Control (available on all i-Limb hands)**

Triggers are specific muscle signals that can be used to access an automated grip. There are 4 potential muscle triggers: hold open, double impulse, triple impulse and co-contraction.

You can use the app to activate and program muscle control.

## **Proximity Control (available only on i-Limb Quantum hands)**

Grip Chips are small Bluetooth devices which activate grips on the i-Limb Quantum hand when the prosthesis is moved close to them or by tapping the Grip Chip.

Before operating a Grip Chip for the first time, pull the plastic tab out of the device to engage battery.

To use Grip Chips, make sure the i-Limb hand is not connected to the My i-Limb app.

- Proximity is accessed by fully opening the i-Limb hand near (15 cm/ 6") a Grip Chip and then relaxing the muscles. Wait until the grip is accessed, which may take up to 3 seconds. To exit a grip accessed via proximity control, give a long open signal.
- Double tapping a Grip Chip activates the grip. Quickly tap a Grip Chip twice, just like double-clicking a computer mouse. The LED on the Grip Chip flashes once when tapping is successful.  
**NOTE:** A 3 second pause between each double tap is required. The pause prevents the grip chip from incorrectly detecting multiple taps in a very short period of time. This can cause the hand to enter a grip and then immediately exit it again.

The i-Limb must be fully open/digits stalled open for proximity or tapping to be successful.

The individual Grip Chips are programmed using the My i-Limb app and can be reprogrammed by the user at any time.

## **Replacing the Grip Chip battery**

The battery is expected to last 6 months. Battery life is reduced if Grip Chip is stored close to the i-Limb hand. Battery is not rechargeable.

- To replace the battery, insert fingernail between upper and lower caps and run around join to release.
- Remove battery. Insert new battery under retaining clip, ensuring writing on battery is facing outwards.
- Replace the caps by pressing upper and lower caps together.

Battery type: CR1616. Replacement batteries are available from customer support upon request.

## I-LIMB HAND COVERS

---

### ***Donning a Cover***

To don a cover, bring i-Limb hand into the position shown in Figure 16 and turn the hand off. The position can also be achieved using the don/doff quick grip.

### ***Donning the i-Limb Skin Active, i-Limb Skin Contour and i-Limb Skin Natural covers:***

1. Line up cover with the fingers of the i-Limb hand, and slide cover downwards.
2. When the fingers are mostly donned, pull the thumb opening over the thumb (**Fig. 17**).
3. Slide the remainder of the cover over the i-Limb hand (**Fig. 18**).
4. Make sure each digit tip is fitted to the cover.
5. Do not pull the cover too tightly over the hand.

### ***Donning the i-Limb Skin Match cover:***

1. Spray the outer surface with isopropyl alcohol (**Fig. 19**).
2. Invert the cover to the level of the finger openings (**Fig. 20**).
3. Ensure the fingers are straight and not bunched.
4. Align the cover finger holes with the digits of the hand (**Fig. 21**).
5. Pull down onto digits of the i-Limb hand.
6. Place thumb opening over thumb digit.
7. Pull the cover carefully over the remainder of the i-Limb hand, without applying too much pressure on the thumb (**Fig. 22**).
8. Check the cover for bunching and make sure the cover tips are fully fitted to the digit tips.
9. Check the function of the i-Limb hand and ensure full opening and closing is possible and the digit tips align.

### ***Doffing all Covers***

1. Bring the i-Limb hand into the same position as for donning and switch off/ to standby.
2. Pull the cover upwards on each digit to release.
3. Carefully ease the full cover off, without applying too much pressure on the thumb.
4. Continue to pull cover upwards, until fully removed.

---

## CLEANING

---

### ***Cleaning instructions***

The i-Limb hand covering, electrode surface and magnetic charger port can be cleaned with a soft damp cloth and mild soap.

Clean the electrode surface regularly.

Clean the cover with isopropyl alcohol once a week to help with disinfection.

**NOTE:** do not submerge the i-Limb hand, cover, electrode surface or magnetic charger port in water. Do not use any strong chemicals.

## TECHNICAL SPECIFICATIONS

<b>Voltage</b>	Range 6 – 8.4 V
<b>Max. Current</b>	5 A
<b>Battery Capacity</b>	Rechargeable lithium polymer 7.4. V (nominal); 2000 mAh capacity; 1,300 mAh capacity
<b>Max hand load limit (static limit)</b>	40 kg/ 88 lb. (extra small) 90 kg/ 198 lb. (small/ medium/ large)
<b>Finger carry load (static limit)</b>	20 kg/ 44 lb. (extra small) i-Limb Quantum: 48 kg/ 106 lb. (small/ medium/ large) i-Limb Ultra and i-Limb Access: 32 kg/ 71 lb. (small/ medium/ large)
<b>Time from full open to full close</b>	0.8 seconds
<b>Expected service life</b>	5 years

Device weight i-Limb Quantum		<b>Extra small</b>	<b>Small</b>	<b>Medium/ large</b>
	QWD	472 g / 1.04 lb.	542 g / 1.19 lb.	558 g / 1.23 lb.
WD		432 g / 0.95 lb.	502 g / 1.1 lb.	518 g / 1.14 lb.
Flexion		572 g / 1.26 lb.	642 g / 1.41 lb.	658 g / 1.45 lb.
Friction		467 g / 1.03 lb.	537 g / 1.18 lb.	553 g / 1.22 lb.

Device weight i-Limb Ultra		<b>Extra small</b>	<b>Small</b>	<b>Medium/ large</b>
	QWD	472 g/ 1.04 lb.	512 g/ 1.13 lb.	528 g/ 1.16 lb.
Note: Titanium digits add an additional 30 g/ 0.07 lbs per hand		432 g/ 0.95 lb.	472 g/ 1.04 lb.	488 g/ 1.08 lb.
Flexion		572 g/ 1.26 lb.	612 g/ 1.35 lb.	628 g/ 1.38 lb.
Friction		467 g/ 1.03 lb.	507 g/ 1.12 lb.	523 g/ 1.15 lb.

Device weight i-Limb Access		<b>Extra small</b>	<b>Small</b>	<b>Medium/ large</b>
	QWD	432 g/ 0.95 lb.	468 g/ 1.03 lb.	478 g/ 1.05 lb.
Note: Titanium digits add an additional 30 g/ 0.07 lb. per hand		392 g/ 0.86 lb.	428 g/ 0.94 lb.	438 g/ 0.97 lb.
Flexion		532 g/ 1.17 lb.	568 g/ 1.25 lb.	578 g/ 1.27 lb.
Friction		427 g/ 0.94 lb.	463 g/ 1.02 lb.	473 g/ 1.04 lb.

## ENVIRONMENTAL CONDITIONS

Do not use, transport or store i-Limb hand outside of the boundaries in the table below:

	<b>Use</b>	<b>Shipping</b>	<b>Extended Storage</b>
<b>Temperature</b>	0 °C to +40 °C	-40 °C to +70 °C	-25 °C to +70 °C
<b>Relative Humidity</b>	10 % to 100 %	10 % to 100 %	10 % to 100 %
<b>Atmospheric pressure</b>	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

## BLUETOOTH MODULES REGULATED INFORMATION

This device contains the following radio frequency transmitters:

Model	Re	Type and Frequency Characteristics	Effective Radiated Power
Bluetooth Low Energy Dual Mode Module Model BR-LE4.0-D2A	FCC  Contains FCC ID: XDULE40-D2  Canada  Contains IC: 8456A-LE4D2  Japan  Contains transmitter with certificate number   [R] 205-160268	(Dual Mode)  Version V2.1 +ED (GFSK + $\pi/4$ DQPSK + 8DPSK) 2402 -2480 MHz  Version V4.0 (GFSK) 2402 -2480 MHz	Adjustable Power (-23 dBm to 10.5 dBm) short to long range

## ELECTROMAGNETIC COMPATIBILITY

**WARNING:** Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

**WARNING:** Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation

In order to regulate the requirements for Electromagnetic Compatibility (EMC) with the aim to prevent unsafe product situations, the BS EN 60601-1-2/ IEC 60601-1-2 standard has been implemented for all Össur Myoelectric Prosthetic Devices. This standard defines levels of electromagnetic emissions for medical devices.

The i-Limb hand is suitable for use in any environment except where immersion in water or any other fluid is possible, or where exposure to highly electrical and/or magnetic fields can occur (e. g. electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

Refer to further guidance below regarding the EMC environment in which the device should be used:

### Guidance and manufacturer's declaration – electromagnetic emissions

Össur Myoelectric Prosthetic Devices is intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that it each are used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Not applicable Battery Powered	Össur Myoelectric Prosthetic Devices use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Össur Myoelectric Prosthetic Devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable Battery Powered	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable Battery Powered	

### Guidance and manufacturer's declaration – electromagnetic immunity

Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.  Portable and mobile RF communications equipment should be used no closer to any part of the Össur Myoelectric Prosthetic Devices, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter
Electrical fast transient/burst IEC 61000-4-4	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered No Cables >3 m
Surge IEC 61000-4-5	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered
Voltage dips, short interruptions and voltage variations on power supply IEC 61000-4-11	Not applicable	Not applicable Battery Powered	Not applicable Battery Powered No Cables >30 m
Power frequency (50/ 60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m 50/60 Hz	Battery Powered

## Guidance and manufacturer's declaration – electromagnetic immunity

Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customers or the users of Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	Not applicable	Not applicable Battery Powered No Cables >3 m	Portable and mobile RF communications equipment should be used no closer to any part of the Össur Myoelectric Prosthetic Devices including cables, than the recommended separation distance calculated from the equation appropriate to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	12 V/m 26 MHz to 1000 MHz  10 V/m 1000 MHz to 2700 MHz	12 V/m 26 MHz to 1000 MHz  10 V/m 1000 MHz to 2700 MHz  1kHz 80% AM	Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ 80 MHz to 800 MHz}$ $d = 2.3 \sqrt{P} \text{ 800 MHz to 2.7 GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters as determined by an electromagnetic site survey <sup>a</sup> should be less than the compliance level in each frequency range <sup>b</sup>  Interference may occur in the vicinity of equipment marked with the following symbol:  

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Össur Myoelectric Prosthetic Devices are used exceeds the applicable RF compliance level above, Össur Myoelectric Prosthetic Devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re- orienting or relocating the Össur Myoelectric Prosthetic Devices

## **Recommended separation distance between portable and mobile RF communications equipment and the Össur Myoelectric Prosthetic Devices**

Össur Myoelectric Prosthetic Devices are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customers or the user of Össur Myoelectric Prosthetic Devices can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Össur Myoelectric Prosthetic Devices as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in Watt	Separation distance according to frequency of transmitter in meters		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## **REPORT A SERIOUS INCIDENT**

### **Important notice to users and/or patients established in Europe:**

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

## DESCRIPTION OF SYMBOLS

BF Applied Part

Manufacturer - YYYY

Manufacturer and year of manufacturing (YYYY)

Consult Instructions for use

Device contains electronic components and/or batteries that should not be disposed of in regular waste

**IP22**

Protected against dust and splashing water from all directions.

**SN**

Serial Number

The unique serial number for i-Limb Quantum devices is a “M” with a 6 digit alpha/numeric number.

The unique serial number for i-Limb Ultra devices is a “U” with a 6 digit alpha/numeric number.

The unique serial number for i-Limb Access is a “A” with a 6 digit alpha/numeric number.

**CE**

European Conformity

**MD**

Medical Device

Class II equipment – provides double isolation to protect against electric shock

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