User manual appendix C

myopro[®]
motion g

Specific instructions for your unique MyoPro 2

myomo[®]
my own motion


PN25623 Rev8

CONTENTS


› MyoPro 2 Motion G Diagram	5
› Finger Support Diagram Close-Up	6
› Putting Your MyoPro 2 On (“Donning”)	7
› Warnings Review	9
› Steps for Putting Your MyoPro 2 On	10
› Taking Your MyoPro 2 Off (“Doffing”)	17
› Warnings Review	19
› Steps for Taking Your MyoPro 2 Off	19
› Using Your MyoPro 2	21
› Control panel	23
› Modes of Use	24
› Start Mode	24
› Elbow Modes	25
› Grasp Modes	27
› Changing Elbow Modes with the Control Panel	29
› Changing Grasp Modes with the Control Panel	29
Dansk	30
Deutsch	60

MYOPRO 2 MOTION G DIAGRAM

Component names and references (not every device includes every component).

FINGER SUPPORT DIAGRAM CLOSE-UP

PUTTING YOUR MYOPRO 2 ON ("DONNING")

PUTTING YOUR MYOPRO 2 ON (“DONNING”)

The following pages include step-by-step instructions for putting the MyoPro 2 on (sometimes referred to as “donning the MyoPro 2”).

WARNINGS REVIEW

BEFORE donning the MyoPro 2, review the warnings below:

-  **WARNING:** DO NOT expose the MyoPro 2 to flame or excessive heat; personal injury may occur.
-  **WARNING:** The MyoPro 2 is a prescription device and should only be used as indicated under clinical supervision of a Provider, or after receiving instruction from a clinical Provider for at-home use.
-  **WARNING:** The MyoPro 2 is only to be used by the person for whom it is prescribed.
-  **WARNING:** Each component of the device is supplied for safety and best performance. Use all provided components for optimal operation.
-  **WARNING:** ONLY use the MyoPro 2 upper limb orthosis on the prescribed elbow, wrist and hand joints. Each MyoPro 2 is built specifically for a left arm or a right arm; never wear the MyoPro 2 on the other arm.
-  **WARNING:** Tight straps may restrict the User’s circulation. Therefore, always check that straps are not too tight throughout the User’s range of motion.
-  **WARNING:** Do not sleep while wearing the device as straps may alter the User’s circulation.
-  **CAUTION:** If you detect fumes, flames, melting of components around the battery or battery charger, or the device is hot to touch, **TURN THE DEVICE OFF IMMEDIATELY** and contact your Provider.

STEPS FOR PUTTING YOUR MYOPRO 2 ON

1. **POWER OFF:** While donning, the MyoPro 2 should be turned OFF; check that it is OFF before starting to put the MyoPor on your arm.
 - ⊙ **NOTE:** Pads and straps should be fully air-dried from the previous use.
 - ▲ **CAUTION:** Do not use a hairdryer to dry components of the MyoPro 2.
2. **PREPARE:** Place the MyoPro 2 in your lap, or on a flat surface that is not slippery.
3. **HAND:** Start with the hand support shell: check that the strap around the hand support shell is loose. Feed your fingers through so the under side of the hand support shell sits in your palm, and your thumb rests in the Thumb Saddle. Be sure the base of your thumb presses against the hard plastic (opposite of the hand strap). Leave the fingers free while you continue to put the rest of the device on – we will secure the finger support as the last step.
 - ⊙ **NOTE:** The motor and battery case should rest on the outside of the arm. If the motor and battery case is aligned along the interior of the arm, you may be donning upside down or on the wrong arm.
4. **WRIST:** Lower your wrist into the multi-articulating wrist (MAW), rotation portion.
5. **FOREARM:** Lower your forearm into the forearm shell.

6. Secure the Forearm Closure.
7. **FINISH LOWER ARM:** Tighten and secure the Velcro hand strap. This will keep the MyoPro 2 in place while you continue to don.

STEPS FOR PUTTING YOUR MYOPRO 2 ON, CONTINUED

8. UPPER ARM:


Wrap the Upper Sensor Cuff around your upper arm, keeping attention on the Biceps and Triceps sensors; when the Upper Sensor Cuff is properly fit, the sensors should be centered over the biceps and triceps muscles. Secure the Upper Cuff Closure and pull to tighten before securing the strap's Velcro.

- ⌚ **NOTE:** the sensor with the **BLUE** Velcro goes on biceps, the sensor with the **RED** Velcro goes on triceps.

Upper Sensor Cuff

Upper Cuff Closure

STEPS FOR PUTTING YOUR MYOPRO 2 ON, CONTINUED

9. FOREARM SENSOR: Slide the Forearm Sensor Cuff onto the forearm. Tuck the Forearm Sensor Cuff's Velcro strap under the forearm bar and secure the Velcro.



- ⓘ **NOTE:** The sensor with the **GREEN** Velcro goes on the inner side of the forearm (aligned with the palm), and the sensor with the **YELLOW** Velcro goes on the outer side of the forearm (aligned with the back of the hand).

STEPS FOR PUTTING YOUR MYOPRO 2 ON, CONTINUED

10. HARNESS: Secure the harness:

- Rest the shoulder pad on the shoulder of your MyoPro 2 arm.
- Wrap the strap under your non-MyoPro 2 arm.
- Secure the Quick Release tab.

STEPS FOR PUTTING YOUR MYOPRO 2 ON, CONTINUED


11. FINGER SUPPORT:

Start with the finger support detached from the Anvil.

Work the Finger Saddle under your fingers while holding the Stem and Loop.

Lift the Stem upward to position the Loop over the Anvil.


You can hang the Loop from the Anvil to take a break or to reposition your other hand for the next step.

STEPS FOR PUTTING YOUR MYOPRO 2 ON, CONTINUED

11. FINGER SUPPORT, continued:


Technique A

Continue to guide the loop over the anvil. The instructions below outline two different techniques for this:

Push down on the end of the Loop to engage the Latch.

Technique b

Lift your fingers from underneath the Finger Saddle.

Since the Stem is rigid, this will cause the end of the loop to rotate over the Latch. (You may also push down on the stem for additional leverage.)

Done!

TAKING YOUR MYOPRO 2 OFF ("DOFFING")

TAKING YOUR MYOPRO 2 OFF (“DOFFING”)

WARNINGS REVIEW

-  **WARNING:** If at any time during the use of this device, you notices any of the following, discontinue use and seek guidance from your Provider.
- Movement does not match the User’s desired motion.
 - Persistent redness, swelling, or skin breakdown (bleeding, chafing, etc.)
 - Rash on the arm, hand, or fingers.
 - Pain associated with wearing the MyoPro 2 orthosis.
 - Unusual noises from the orthosis (popping, clicking, etc.)
 - Smells from the orthosis (smoking, burning plastic, etc.)
 - Odor from the orthosis (sour smells or other indications of bio-contamination.)

STEPS FOR TAKING YOUR MYOPRO 2 OFF


1. Turn the MyoPro 2 OFF.
2. Place the MyoPro 2 in your lap, or on a flat surface that is not slippery.
3. Undo the harness Quick Release tab.
4. Release the Finger Support (see Page 17 for images).
5. Loosen the hand strap.
6. Remove the Forearm Sensor Cuff from your forearm.
7. Undo the Upper Cuff Closure.
8. Undo the Forearm Closure.
9. Gently slide your hand back through the loosened hand strap and Hand Support Shell
10. Rest the MyoPro 2 in a clean, dry location to ensure it will fully air-dry before next use.
11. Place the MyoPro 2 into the carrying bag for travel or storage.

-  **NOTE:** The MyoPro 2 should remain OFF any time the MyoPro 2 is not being worn.


STEPS FOR TAKING YOUR MYOPRO 2 OFF, CONTINUED

Push/pull the slider forward to release the loop.

The loop will slide over the Anvil as soon as it is free.


Remove the "W", Stem and Loop from under your fingers and store in the MyoPro 2's carrying bag.

USING YOUR MYOPRO 2

CONTROL PANEL

See the *User Manual – Sections for All MyoPro 2 Devices* document for a comprehensive Control Panel Buttons & Lights Diagram.

MODES OF USE

 **WARNING:** Use caution when using the device in shoulder positions where it is possible for the User to hit him or herself.

The MyoPro 2 may be used in a variety of **modes**. The different modes allow you to operate the MyoPro 2 by tensing or relaxing different muscles. You may find that you want to use different muscles – and therefore different modes – when performing specific tasks with your MyoPro 2. Modes may be changed using the buttons on the control panel.

When the device is fit to you, your Provider will optimize the sensitivity settings within each mode depending on the strength of your EMG signal at that time. If over time you feel that the MyoPro 2 is not assisting your arm or hand adequately, please contact your Provider; the sensitivity settings may need to be readjusted.

START MODE

When you first power on the MyoPro 2, the elbow will be in Standby Mode. This setting may be changed by your Provider.

ELBOW MODES

The elbow has four distinct modes of use when the MyoPro 2 is powered on.

1. Standby mode
2. Biceps mode
3. Triceps mode
4. Dual mode

These are explained in the following pages.

GRASP MODES

The grasp has four distinct modes of use when the MyoPro 2 is powered on.

1. Standby mode
2. Open mode
3. Close mode
4. Dual mode

These are explained in the following pages.

 **NOTE:** You may experience a slight vibration of the finger support when the grasp reaches a fully closed or fully open position. This would be caused by the hand motor turning on and off as it detects the limits of its range of motion. This is normal operation, and nothing to worry about.

ELBOW MODES

<p>Standby mode</p>	<p>In this mode, neither the Biceps or Triceps Light will be illuminated. The elbow motor will not respond to your EMG signal from either muscle group.</p> <p>Though the MyoPro 2 is powered on, no assistance is being given to your elbow, and the arm will neither flex nor extend. The sensors will be reading your EMG signal from both biceps and triceps, but the motor response is paused.</p> <p>Actions:</p> <ul style="list-style-type: none"> • This can be useful as a resting mode if straps or sensors need to be adjusted, or if you want to take a short break without turning the device off. 	
<p>Biceps mode</p>	<p>In this mode, the Biceps Light will be illuminated. The elbow motor will respond to your biceps EMG signal.</p> <p>Actions:</p> <ul style="list-style-type: none"> • When you relaxes your biceps, the device will extend. • When you contracts your biceps, the device will flex. 	

ELBOW MODES, CONTINUED

<p>Triceps mode</p>	<p>In this mode, the Triceps Light will be illuminated. The elbow motor will respond to your triceps EMG signal.</p> <p>Actions:</p> <ul style="list-style-type: none"> • When you relax your triceps, the elbow will flex. • When you contracts your triceps, the elbow will extend. 	
<p>Dual mode</p>	<p>In this mode, both the Biceps and Triceps Light will illuminated. The elbow motor will respond to your biceps <i>and</i> triceps EMG signal.</p> <p>The elbow will only respond to muscle contraction (not relaxation, as in other modes) to assist you with active flexion and active extension of your elbow.</p> <p>You must engage both muscle groups in order for your arm to move in both directions. The MyoPro 2 will respond to whichever muscle signal is stronger at a given moment in time. In order for Dual mode to be most effective, you should be able to isolate the desired muscle (biceps or triceps) and have a limited number of co-contractions prior to using Dual mode. When used well, this mode should allow you greater ability to hold an arm position at a mid-point between fully flexed and fully extended.</p> <p>Actions:</p> <ul style="list-style-type: none"> • When you contract your biceps, the elbow will flex, as long as your biceps signal is greater than your triceps signal. • When you contract your triceps, the elbow will extend, as long as your triceps signal is greater than your biceps signal. 	

GRASP MODES


<p>Standby mode</p>	<p>In this mode, neither the Close or Open Light will be illuminated. The grasp motor will not respond to your EMG signal.</p> <p>Though the MyoPro 2 is powered on, no assistance is being given to your hand, and the grasp will neither open nor close. The sensors will be reading your EMG signals from both wrist flexor and extensor muscle groups, but the motor response is paused.</p> <p>Actions:</p> <ul style="list-style-type: none"> • This can be useful as a resting mode if straps or sensors need to be adjusted, or if you want to take a short break without turning the device off. 	
<p>Close mode</p>	<p>In this mode, the Close Light will be illuminated. The grasp motor will respond to your wrist flexor EMG signal.</p> <p>Actions:</p> <ul style="list-style-type: none"> • When you relax your wrist flexors, the grasp will open. • When you contract your wrist extensors, the grasp will close. 	

GRASP MODES, CONTINUED

<p>Open mode</p>	<p>In this mode, the Open Light will be illuminated. The grasp motor will respond to your wrist extensor EMG signal.</p> <p>Actions:</p> <ul style="list-style-type: none"> • When you relax your wrist extensors, the grasp will close. • When you contract your wrist extensors, the grasp will open. 	
<p>Dual mode</p>	<p>In this mode, both the Close and Open Light will illuminated. The grasp motor will respond to your wrist flexor <i>and</i> wrist extensor EMG signal.</p> <p>The grasp will only respond to muscle contraction (not relaxation, as in other modes) to assist you with active close and active open of your hand.</p> <p>You must engage both muscle groups in order for your hand to move in both directions. The MyoPro 2 will respond to whichever muscle signal is stronger at a given moment in time. In order for Dual mode to be most effective, you should be able to isolate the desired muscle group (wrist flexors or wrist extensors) and have a limited number of co-contractions prior to using Dual mode. When used well, this mode should allow you greater ability to hold a hand position at a mid-point between fully closed and fully open.</p> <p>Actions:</p> <ul style="list-style-type: none"> • When you contract your wrist flexors, the grasp will close, as long as your wrist flexor signal is greater than your wrist extensor signal. • When you contract your wrist extensors, the grasp will open, as long as your wrist extensor signal is greater than your wrist flexor signal. 	

CHANGING ELBOW MODE WITH THE CONTROL PANEL

The elbow mode button cycles through the 4 elbow modes, moving to the next mode each time you press the mode button.

When the MyoPro 2 is turned ON, the elbow will be in Standby mode unless it has been changed by your Provider.	
	Press the elbow mode button one time to put the elbow into Biceps mode
	Press the elbow mode button again (a second time) to put the elbow into Triceps mode
	Press the elbow mode button again (a third time) to put the elbow in Dual mode
	Press the elbow mode button again (a fourth time) to return the elbow Standby mode

Repeat the above cycle to scroll through and change to another mode at any point during your MyoPro 2 use.

CHANGING GRASP MODE WITH THE CONTROL PANEL

The grasp mode button cycles through the 4 grasp modes, moving to the next mode each time you press the mode button.

When the MyoPro 2 is turned ON, the grasp will be in Standby mode unless it has been changed by your Provider.	
	Press the grasp mode button one time to put the grasp into Close mode
	Press the grasp mode button again (a second time) to put the grasp into Open mode
	Press the grasp mode button again (a third time) to put the grasp in Dual mode
	Press the grasp mode button again (a fourth time) to return the grasp Standby mode

Repeat the above cycle to scroll through and change to another mode at any point during your MyoPro 2 use.

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